

VasoStat™ HEMOSTASIS DEVICE INSTRUCTIONS FOR USE

CAUTION: Rx only: U.S. federal law restricts this device to sale to or on the order of a licensed health care practitioner.

DEVICE DESCRIPTION: The VasoStat™ Hemostasis Device is used to augment hemostasis through a ratcheting pressure mechanism.

INTENDED USE: The VasoStat™ Hemostasis Device is indicated for use by medical professionals to promote hemostasis following a catheterization or other puncture into a blood vessel in a patient's arm or lower leg, including radial artery catheterization, pedal or tibial artery catheterization, arterial or venous line removal, hemodialysis, and in patients on anticoagulation therapy.

CONTRAINDICATIONS: Allergy to medical adhesives.

PRECAUTIONS:

- The device is intended for single use only. Do not sterilize and/or reuse this device.
- Do not use this product if there is doubt as to whether the product is sterile.
- Store in a dark, dry, cool place. Avoid extended exposure to light.
- Upon removal of the product from packaging, inspect product to ensure no damage has occurred.

POTENTIAL ADVERSE EVENTS:

- Undesired bleeding at the puncture site could occur if care is not used when securing the device at puncture site
- Hematoma
- Radial artery or tibial artery thrombosis if the device is tightened excessively without checking for the presence of a pulse and/or hand/foot perfusion distal to the site of device application
- Vascular thrombosis if the device is tightened excessively over a hemodialysis access without checking for the presence of a thrill
- As in common with skin adhesives, slight irritation of the skin may occur with prolonged use

INSTRUCTIONS FOR USE:


1. Remove paper strips on the base of the VasoStat™ Hemostasis Device to expose the adhesive surfaces. Ensure plunger of device is in the 'up' position. If needed, gently squeeze the wings of the plunger together to move plunger to the 'up' position.

2. Align central plunger of device above vascular puncture site while vascular sheath or catheter is still in place. If needed, partially withdraw the vascular sheath or catheter to provide sufficient room for positioning and securement of the device.

3. Secure adhesive wings of VasoStat™ device to patient's skin with manual pressure. Ensure the adhesive pads completely contact the patient's skin. Next, remove the paper strips and apply the SeeCure™ adhesive pad over the wings of the VasoStat™ device and surrounding skin.

4. While stabilizing the base of the device with index and middle fingers, depress the plunger of VasoStat™ device with thumb pressure until it is just contacting the patient's skin.

5. Remove vascular sheath or catheter while simultaneously depressing plunger of VasoStat™ to apply desired amount of compression force to the puncture site.




6. If bleeding is not controlled, apply additional pressure by further depressing the plunger against the puncture site.

7. Ensure a pulse can still be felt beyond the device so that an excessive amount of force is not being applied. If unsure, loosen the device (see steps 8-9 below) and check for the presence of a pulse before re-tightening the plunger.

8. After sufficient time to achieve hemostasis, release the plunger by squeezing the wings on each side to disengage them from the locking mechanism within the device.

9. Withdraw the central plunger until it no longer contacts the skin and/or no longer produces pressure on the skin. Verify hemostasis has occurred by visually inspecting the puncture site underneath the plunger.

10. Remove the VasoStat™ Hemostasis Device and SeeCure™ pad by gently peeling the adhesive surfaces free from the patient's skin. Isopropyl alcohol swab(s) may also be used to assist in removal. Discard the device in appropriate receptacle. A sterile dressing may be applied to the puncture site.



HOW SUPPLIED:

The VasoStat™ Hemostasis Device and SeeCure™ adhesive pad are supplied sterile in peel open packages, intended for one-time use. Sterile if package is unopened or undamaged.

Please report any defects or adverse events to 877-466-0109.

 **Forge Medical, Inc.**
 2436 Emrick Boulevard
 Bethlehem, Pennsylvania, US, 18020
 + 1 877 466 0109

EC REP

MedQ Consultants B.V.
 NL Kloosterweg 1
 6412 CN Heerlen
 + 31 45 303 0006



Do not resterilize



Non Latex



Do not reuse



Do not use if package or seal is damaged



Consult Instructions for Use



Contents: 1 Device



Manufacturer



Use by

STERILE R

Sterilized by gamma radiation



CE Marking



Lot number

EC REP

Authorized Representative in the European Community



Telephone number



Facsimile number

